

Product datasheet

Anti-GKAP1 antibody - C-terminal ab179412

1 Image

Overview

| | |
|----------------------------|--|
| Product name | Anti-GKAP1 antibody - C-terminal |
| Description | Rabbit polyclonal to GKAP1 - C-terminal |
| Host species | Rabbit |
| Tested applications | Suitable for: WB |
| Species reactivity | Reacts with: Human Predicted to work with: Rat, Rabbit, Horse, Guinea pig, Cow, Cat, Dog |
| Immunogen | Synthetic peptide within Human GKAP1 aa 317-366 (C terminal). The exact sequence is proprietary. (NP_079487). Sequence: SIKNELTIQVTSLHAALEQERSKVKVLQAELAKYQGGRKG KRNSEDQCR Database link: Q5VSY0 Run BLAST with Run BLAST with |
| Positive control | HCT15 whole cell lysate. |

Properties

| | |
|-----------------------------|---|
| Form | Liquid |
| Storage instructions | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle. |
| Storage buffer | Constituents: 98% PBS, 2% Sucrose |
| Purity | Immunogen affinity purified |
| Clonality | Polyclonal |
| Isotype | IgG |

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab179412 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


| Application | Abreviews | Notes |
|-------------|-----------|---|
| WB | | Use a concentration of 1 µg/ml. Predicted molecular weight: 42 kDa. |

Target

Sequence similarities Belongs to the GKAP1 family.

Cellular localization Golgi apparatus.

Images



Anti-GKAP1 antibody - C-terminal (ab179412) at 1 µg/ml + HCT15 whole cell lysate at 10 µg

Predicted band size: 42 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors