

Product datasheet

Anti-DALRD3 antibody ab201209

2 Images

Overview

<b>Product name</b>	Anti-DALRD3 antibody
<b>Description</b>	Rabbit polyclonal to DALRD3
<b>Host species</b>	Rabbit
<b>Tested applications</b>	<b>Suitable for:</b> WB, IHC-P
<b>Species reactivity</b>	<b>Reacts with:</b> Human
<b>Immunogen</b>	<p>Recombinant fragment corresponding to Human DALRD3 aa 194-396. (BC054493). Sequence:</p> <p>HALEELTSANDGRTLSPGILGRLCLKELVEEQGRTAGYDP            NLDNCLVTED            LLSVLAELQEALWHWPEDSHPLAGASDTGTGGCLVVH            VVSCEEEFQQQK            LDLLWRKLVDKAPLRQKHLICGPVKVAGAPGLMTAPEYY            EFRHTQVCKA            SALKHGGDLAQDPAWTEIFGVLSVATIKFEMLSTAPQSQL            FLALADSSIS TKG</p> <p>Database link: <a href="#">Q5D0E6-2</a>  <a href="#">Run BLAST with</a>  <a href="#">Run BLAST with</a></p>
<b>Positive control</b>	HepG2 cell lysate; Human fetal colon tissue.

Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	<p>pH: 7.2            Preservative: 0.02% Sodium azide            Constituents: 1% BSA, 98% PBS</p>
<b>Purity</b>	Immunogen affinity purified
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG


## Applications

**The Abpromise guarantee** Our [Abpromise guarantee](#) covers the use of ab201209 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/200 - 1/1000. Predicted molecular weight: 55 kDa.
IHC-P		1/100 - 1/500.


## Target

## Images



Anti-DALRD3 antibody (ab201209) at 1/1000 dilution + HepG2 cell lysate

**Predicted band size:** 55 kDa



Immunohistochemical analysis of formalin-fixed, paraffin-embedded Human fetal colon tissue labeling DALRD3 with ab201209 at 1/100 dilution.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-DALRD3 antibody (ab201209)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors