

Product datasheet

Anti-12L HETE antibody ab2275

Overview

Product name	Anti-12L HETE antibody
Description	Rabbit polyclonal to 12L HETE
Host species	Rabbit
Specificity	This antibody cross reacts with Methyl ester-12-HETE (20%), 8-HETE (8.7%) and 12(L)-HHT (8.6%).
Tested applications	Suitable for: RIA, ELISA
Species reactivity	Reacts with: Species independent
Immunogen	Pure 12(L)-HETE conjugated to human serum albumin
General notes	Principal Name - 12(L)- Hydroxeicosantetraenoic Acid.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	unknown
Purity	Whole antiserum
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab2275 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
RIA		Use at an assay dependent concentration.
ELISA		Use at an assay dependent concentration.

Target

Relevance

12R-hydroxyeicosatetraenoic acid (12R-HETE) is produced following the catalytic conversion of arachidonic acid by 12R-lipoxygenase. Boeglin et al. (1998) provided mechanistic evidence for a lipoxygenase route to 12R-HETE in human psoriatic tissue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors