

Product datasheet

Anti-2,3 dinor 6 keto Prostaglandin F1 alpha antibody ab2276

1 References

Overview

Product name	Anti-2,3 dinor 6 keto Prostaglandin F1 alpha antibody
Description	Sheep polyclonal to 2,3 dinor 6 keto Prostaglandin F1 alpha
Host species	Sheep
Specificity	This antibody cross reacts with 6-keto-PGF1 alpha (89.14%), PGF1 alpha (13.21%), PGF2 alpha (4.99%) and 13,14-dihydro-PGF1 alpha (2.55%). The antibody will detect 2,3 dinor-6-keto-Prostaglandin F1 alpha in any species.
Tested applications	Suitable for: ELISA, RIA
Immunogen	Highly pure 2,3-dinor-6-keto-PGF1 alpha conjugated to BSA

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	unknown
Purity	Whole antiserum
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab2276 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		Use at an assay dependent concentration.
RIA		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors