

Product datasheet

Anti-Prostaglandin F2 alpha antibody ab2396

Overview

Product name	Anti-Prostaglandin F2 alpha antibody
Description	Rabbit polyclonal to Prostaglandin F2 alpha
Host species	Rabbit
Specificity	This antibody cross-reacts with Prostaglandin F1 alpha (100%), Prostaglandin E1 (1.1%) and 6-keto-Prostaglandin F1 alpha (1.1%).
Tested applications	Suitable for: ELISA
Immunogen	Prostaglandin F2 alpha conjugated to a carrier molecule.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	unknown
Purity	Whole antiserum
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab2396 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		

Application notes ELISA: 1/5000.

Not yet tested in other applications.
Optimal dilutions/concentrations should be determined by the end user.

Target

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors