

Product datasheet

Anti-THTPA antibody ab57728

[2 Images](#)

Overview

Product name	Anti-THTPA antibody
Description	Mouse monoclonal to THTPA
Host species	Mouse
Tested applications	Suitable for: WB, IHC-P
Species reactivity	Reacts with: Human
Immunogen	Recombinant full length protein, corresponding to amino acids 1-231 of Human THTPA

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: None PBS, pH 7.2
Purity	Protein G purified
Clonality	Monoclonal
Isotype	IgG1
Light chain type	kappa

Applications


The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab57728 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 25 kDa.
IHC-P		Use a concentration of 3 µg/ml.

Target


Relevance	Thiamine triphosphatase (THTPA) catalyzes the hydrolysis of thiamine triphosphate (THTP) to thiamine diphosphate (THDP), which is the major cellular form of thiamine (vitamin B1).
Cellular localization	Cytoplasm.

Images



THTPA antibody (ab57728) used in immunohistochemistry at 3ug/ml on formalin fixed and paraffin embedded human small Intestine.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-THTPA antibody (ab57728)



THTPA antibody (ab57728) at 1ug/lane + A-431 cell lysate at 25ug/lane.

Western blot - Anti-THTPA antibody (ab57728)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors