



Product datasheet

Anti-GSK3 beta (phospho S9) antibody ab9769

[4 References](#) [2 Images](#)

Overview

Product name	Anti-GSK3 beta (phospho S9) antibody
Description	Rabbit polyclonal to GSK3 beta (phospho S9)
Host species	Rabbit
Specificity	Gives strong signal to GSK 3 beta in westerns at the estimated mol weight (49 kD). May show very weak reactivity to GSK 3 alpha.
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human
Immunogen	Synthetic peptide: RRTT _p SFAE, corresponding to amino acids 4-12 of Human GSK3 beta. <div>  Run BLAST with  Run BLAST with </div>

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, 0.01% Sodium Azide, pH 7.2
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab9769 in the following tested applications.
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/1000.

Target

Function

Participates in the Wnt signaling pathway. Implicated in the hormonal control of several regulatory proteins including glycogen synthase, MYB and the transcription factor JUN. Phosphorylates JUN at sites proximal to its DNA-binding domain, thereby reducing its affinity for DNA. Phosphorylates MUC1 in breast cancer cells, and decreases the interaction of MUC1 with CTNNB1/beta-catenin. Phosphorylates CTNNB1/beta-catenin. Phosphorylates SNAI1. Plays an important role in ERBB2-dependent stabilization of microtubules at the cell cortex. Prevents the phosphorylation of APC and CLASP2, allowing its association with the cell membrane. In turn, membrane-bound APC allows the localization of MACF1 to the cell membrane, which is required for microtubule capture and stabilization. Phosphorylates MACF1 and this phosphorylation inhibits the binding of MACF1 to microtubules which is critical for its role in bulge stem cell migration and skin wound repair.

Tissue specificity

Expressed in testis, thymus, prostate and ovary and weakly expressed in lung, brain and kidney.

Sequence similarities

Belongs to the protein kinase superfamily. CMGC Ser/Thr protein kinase family. GSK-3 subfamily. Contains 1 protein kinase domain.


Post-translational modifications

Phosphorylated by AKT1 and ILK1. Activated by phosphorylation at Tyr-216.

Cellular localization

Cytoplasm. Nucleus. Cell membrane. The phosphorylated form shows localization to cytoplasm and cell membrane. The MEMO1-RHOA-DIAPH1 signaling pathway controls localization of the phosphorylated form to the cell membrane.

Images



Rabbit polyclonal to phospho GSK beta (Ser 9) at 1/1000 on 293T cell extract shows a band at 49 kD.

Cells were serum starved for 24 hours:

Lane 1 : control

Lane 2 : treated with IGF-1 (100 ng/ml) for 20 min

Lane 3 : pre-treated with 10 uM LY294002 (selective PI3K inhibitor) and treated with IGF-1 (100 ng/ml) for 20 min


Rabbit polyclonal to phospho GSK beta (Ser 9) at 1/1000 on 293T cell extract shows a band at 49 kD.

Cells were serum starved for 24 hours:

Lane 1 : control

Lane 2 : treated with IGF-1 (100 ng/ml) for 20 min

Lane 3 : pre-treated with 10 uM LY294002 (selective PI3K inhibitor) and treated with IGF-1 (100 ng/ml) for 20 min



IHC-P of ab9769 at 1:100 dilution detecting phosphorylated Huntingtin in human brain cerebellum. Positive cytoplasmic staining is observed in neurons.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-GSK3 beta (phospho S9) antibody (ab9769)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors